

GO GO[®]

ELITE TRAVELLER GEN 3

PLATINUM

FDA Class II Medical Device

Featuring a new deck design for superior comfort and **increased leg room!**



New and improved deck design for increased comfort and leg room.



Intuitive control panel framed by a soft-touch delta tiller handle.



EZ Turn technology provides unrivaled tight turning without sacrificing stability.

TRAVEL THE WORLD IN COMFORT AND STYLE

The new **Go Go Elite Traveller[®] Gen 3 Platinum** takes comfort and convenience to the next level with increased leg room and exclusive luxury features. Style upgrades include color accents on the front shroud and wheels, while the front basket has been upgraded to a larger size for increased storage of your personal belongings.

You'll ride with pride, knowing your scooter was built in the USA from foreign and domestic parts. Ask your local authorized **Pride[®] Dealer** for a test drive today and experience the **platinum** difference.



SC443EMS

Go Go Elite Traveller® Gen 3 Platinum

FDA Class II
Medical Device



INCLUDED COLOR PANELS (MATTE)



Cardinal Red

ACCESSORY COLOR PANELS (MATTE)



Black



Iceberg
Blue



Robin's
Egg Blue



Sugar
Plum



Chell-o

HCPCS Code: K0899

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 3 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 5 Battery weight may vary based on manufacturer.
- 6 Dimensions listed are for scooter only and do not include seats or arms. Overall width may vary based on seat options and accessories.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

AUTHORIZED DEALER

FEATURES

- **New sloped deck** design for improved comfort and leg room
- **Front and rear pathway lights**
- **EZ Turn** technology provides unrivaled tight turning without sacrificing stability
- **Feather-touch disassembly** allows for easy separation into 5 lightweight pieces for transport or storage
- Exclusive **Charge360** magnetic charging port requires minimal dexterity and lights up to signal a solid connection
- Tiller angle-adjusting foot pedal
- Interchangeable shroud color panels (comes standard with Cardinal Red)
- Large front basket (standard)



ACCESSORIES *

- Additional shroud panels
- Cup holder
- Military patches
- Oxygen holder
- Rear basket
- Rear mount 3" double cane/crutch holder
- Rear view mirror
- Saddle bag
- Safety flag
- Single rear mount 3" cane/crutch holder
- Under seat storage
- Walker holder
- Weather cover

* Only one rear accessory can be used at a time.

SPECIFICATIONS

WEIGHT CAPACITY	300 lbs.	
MAXIMUM SPEED²	Up to 4.23 mph	
DIMENSIONS^{4,6}	Width: 19.45"	Length: 41.5"
TURNING RADIUS⁴	39.2"	
GROUND CLEARANCE	1.29" at motor	3.5" under deck
WEIGHT⁵	83.6 lbs. without batteries 103 lbs. with 12 Ah batteries 111 lbs. with 18 Ah batteries	
RANGE^{2,3}	12 Ah: Up to 7.1 miles at 300 lbs. Up to 8.9 miles at 220 lbs. 18 Ah: Up to 11.3 miles at 300 lbs. Up to 14.1 miles at 220 lbs.	

WARRANTY

FRAME: Lifetime limited
ELECTRONICS: 2-year limited

DRIVETRAIN: 2-year limited
BATTERY: 6 months



pridemobility.com

GO GO
TRAVEL MOBILITY®